Operating Instructions 使用说明书

High-speed Micro Centrifuge 高速微量离心机

CF221

Contents

1. Declaration of Conformity	02
2. Safety Instructions	02
3. Operating Instructions	06
4. Unpacking	10
5. Product Introduction	10
6. Specifications	13
7. Installation Requirements	13
8. Instructions for the Rotor and Tubes	14
9. Operation Settings	16
10. Troubleshooting	20
11. Cleaning and Maintenance	21
12. Warranty	23

Please read the operating instructions in full before starting up and follow the safety instructions.
The appearance and specifications are subject to change without notice.

1. Declaration of Conformity

We declare under our sole responsibility that this product is in compliance with the regulations:

2014/35/EU, 2006/42/EC, 2014/30/EU, 2011/65/EU, EN 61010-1, EN 61010-2-010, EN 61010-2-051, EN 61326-1, EN 60529, EN ISO 12100.

2. Safety Instructions

General Information

- Please read the operating instructions completely before using the device.
- Keep the operating instructions in a place where it can be accessed by users.
- Ensure that only trained staff can work with the device.
- Follow the safety instructions, guidelines, occupational health, safety and accident prevention regulations.

Explanation of Symbols

This symbol indicates prohibited operational actions.

This symbol indicates important mandatory actions to follow.

- ❗ • Prior to usage, please strictly follow the operating steps mentioned in this manual. Failure to do so may damage the device.
- Only trained staff can operate or repair this device.
- Position the device on a secure and steady work surface to prevent any potential harm to the device or personal injury.
- Prior to usage, ensure that the power supply voltage matches the specified requirement indicated on the device nameplate. Neglecting to do so could lead to fire or electric shock.
- The device uses a three-pin plug and the third pin is a ground pin. Don't remove the grounding prong of the power cord and ensure that the power plug is securely connected to the socket to avoid the risk of fire.
- In the event of a gas leak, such as inflammable or explosive gas, close the gas valve, open doors and windows for ventilation, and do not plug or unplug electrical appliances, as this could potentially cause explosions and fires.
- When unplugging the centrifuge plug from the power socket, grasp the plug head instead of pulling the cord. Pulling the cord directly with your hand may result in electric shock or a short circuit, leading to a fire.
- If the device is not working properly, please disconnect the power. Operating the device under abnormal conditions may cause electric shock, leading to a fire.

- Ensure the device is disconnected to the power supply before starting repairs or maintenance to prevent electric shock or personal injury.
- Ensure that no suspended particles from the centrifuge and its surroundings are inhaled during maintenance and cleaning, as this could potentially pose health hazards.
- Disconnect the power supply if the device won't be used for an extended period of time to prevent the power cord from aging.
- The disposal of scrapped centrifuges should be carried out by certified personnel.
- When restarting the centrifuge after a power outage, it is important to check the settings of the centrifuge. Changes in settings could potentially affect the results.
- Be careful when transporting the device to avoid any personal injury or causing damage to the device.
- The device should be used in the open environment of a well-ventilated laboratory or a fume hood.

- Don't place the centrifuge in damp or near water sources, as it may result in electric shock.
- Don't pour water directly onto the centrifuge, as it may cause electric shock or short circuits.
- Don't touch any electrical components of the centrifuge, such as power plug, or any switches with wet hands, as it may cause electric shock.

-  Users are prohibited from disassembling, repairing, or modifying the centrifuge themselves, as improper operation may lead to fires or injuries.
- Don't place flammable, explosive, or volatile substances in the centrifuge, and refrain from using flammable aerosols near the centrifuge to avoid explosions or fires.
- Don't insert metal objects such as nails or wires into any holes or gaps of the centrifuge or into any exhaust vents used for internal air circulation, as contact between these objects and moving parts may cause electric shock or injury.

Disclaimer

The following situations may affect the normal use of the equipment. We shall not responsible for any property loss or personal injury caused by below situations:

- Don't follow the operating instructions
- Use of the equipment for purposes other than specified in this manual
- Use of accessories and consumables not recommended by the manufacturer
- Maintenance or repair of the equipment by uncertified personnel
- Unauthorized modification of the equipment by the user

3. Operating Instructions

Danger! Explosion hazard!

- Don't operate this equipment indoors where there is a risk of explosion.
- Don't use this equipment to centrifuge explosive or highly reactive substances.
- Don't use this equipment to process materials that produce explosive gases.

Warning!

Electric shock due to damaged equipment or power cord.

- Before starting the equipment, ensure that both the equipment and power cord are undamaged.
- Only start the equipment that has been correctly installed or repaired.
- In hazardous situations, disconnect the equipment from the power supply. Remove the power plug from the equipment or power socket.

Warning!

Hazards due to damaged accessories!

- Ensure all accessories are in good condition.
- Check if the accessories are damaged before use. If yes, replace the accessories.
- Don't use the accessories that exceed the service life.

Caution! Hazards due to incorrect accessories!

- Using accessories that aren't recommended by the manufacturer may affect the safety, normal function and accuracy of the equipment. Damages caused by improper accessories is not covered by the warranty.

Device damage caused by contacting liquid

- Turn off the device.
- Disconnect the power supply.
- Carefully clean the equipment and its components according to the cleaning and maintenance instructions mentioned in this manual.

Condensation can lead to damage of electrical components

- When transporting the equipment from a cooler environment to a warmer one, condensation may form on the chamber wall. Please wait for at least 4 hours before reboot the equipment.
- The equipment should be placed in a location where the power supply can be easily disconnected.
- Colliding with or moving the device while it is running may cause damage. Severe damage to the centrifuge or rotor may occur if the rotor collides with the chamber wall.

Rotor and rotor lid are not secured properly

- Tighten the rotor and the rotor lid before running.
- Improper securing of the rotor and rotor lid may cause abnormal noise in the centrifuge. Press this  button immediately to stop the device.

The tubes are asymmetrically loaded

- To ensure symmetric loading, tubes that are arranged opposite each other must be of the same type and contain the same sample quantity.
- To make sure that tubes are symmetrically loaded, we recommend using electronic scale to measure the filled sample tubes.

Rotor Overload

- When the device reaches the maximum speed and maximum loading capacity, the density of the sample tubes must not exceed 1.2g/ml.

The rotor lid is corroded by chemical substances.

- If the rotor lid comes into contact with organic solvents, clean it immediately.
- Regularly inspect the rotor lid for damages or cracks.
- If there are cracks or it turns milky white, replace it immediately.

The rotor is damaged by corrosive chemicals.

- The stability of the rotor may be affected by corrosive chemicals.

- Avoid contacting corrosive chemicals, including strong alkalis, strong acids, solutions containing mercury, copper, and other heavy metal ions, organic halides, concentrated hydrochloric acid solutions, and phenols.
- If the rotor is contaminated with corrosive chemicals, clean it immediately with a neutral detergent and pay attention to clean the grooves of the rotor.

Rotor falls off due to improper operation.

- Always use both hands to carry the rotor when removing or placing it to prevent it from falling off.

Damaged tubes.

- Don't use damaged tubes for centrifugation as it may damage to the device and the sample. Visually inspect all centrifuge tubes for damage before use.

The tube lid is open.

- Unclosed tube lids may break during centrifugation, causing damage to the rotor and centrifuge. Please make sure that all tube lids are closed before running.

Organic solvents can damage plastic centrifuge tubes.

- Using organic solvents (such as phenol) can weaken the strength of plastic centrifuge tubes, potentially causing damage to the tubes.

4. Unpacking

Unpack the package carefully and examine the contents for potential damages that may have occurred during transit. If you observe discernible damage, please promptly contact your local seller for assistance. In the event of any discernible damage to the device, refrain from connecting the device to a power source to avoid any potential risks or further harm.

Packing List:

Name	Quantity
Centrifuge	1
Power cord	1
Emergency lid release pin	1
Wrench	1
User manual	1
Warranty card	1
Certificate	1

5. Product Introduction

Product overview


CF221 High-speed Micro Centrifuge is a versatile laboratory instrument designed for rapid separation and extraction of components from mixtures of liquids and solids. By utilizing centrifugal force, it effectively separates particles based on their density, enabling efficient sample processing in various fields, including life sciences, biology, medical research, pharmaceuticals, geology, agriculture, forestry, petrochemicals, light industry, and

metal particles analysis, etc.

Features

- Maximum centrifugal force (RCF) reaches up to 21,200xg
- Large TFT touchscreen with an intuitive user interface for easy operation
- Motor-driven lock lid design enhances safety, enabling easy lid closure with gentle pressure
- Silent operation promotes a safe, stress-free environment
- Multiple protection function, over speed protection, system protection and lid locking protection
- Remarkable versatility with 4 different rotors: 24x1.5/2.0mL tubes, 12x5mL tubes, 8x4x0.2mL PCR strips/tubes (for 4x0.2mL PCR strips or 32x0.2mL PCR tubes)
- The rotor chamber is made of fully enclosed steel structure, ensuring safe centrifugation

Product structure

①	Lid	⑥	Emergency lid release hole
②	Rotor chamber	⑦	View port
③	Sealing ring	⑧	Fan
④	Operation panel and display	⑨	Power switch
⑤	Indicator	⑩	USB interface

6.Specifications

Model	CF221
Max. Speed (rpm)	15,000
Max. RCF (xg)	21,200
Speed accuracy	±10rpm
Time setting range	1s~99min59sec or continuous run function
Noise	< 60dB(A)
Power supply	AC110V/220V 50-60Hz
Wattage	230W
Weight	17.5kg
Dimension (L*W*H)	382*266*256.8mm

Note: The appearance and technical parameters of the device are subject to change without notice.

Model No.	Rotor Capacity
CF221-01 (Standard)	24x1.5 / 2.0ml tubes
CF221-03 (Optional)	12x5ml tubes
CF221-04 (Optional)	4x8x0.2ml PCR strips/tubes

7. Installation Requirements

Installation Environment

To ensure that the device can meet the normal operation and optimal performance, please abide by the following requirements:


- Indoor use only
- Ambient temperature: 5~35°C
- Avoid strong vibrations, high temperatures, and direct sunlight in the vicinity.

- Avoid mechanical shaking or vibration.
- Position the device on an secure and stable work platform in an area free from corrosive gases, smoke, direct sunlight, excessive air currents, or strong magnetic field disruptions.
- After installation, the four rubber feet of the device should bear even force.
- The distance between the back and side of the device and the wall should not be less than 30cm to ensure well ventilation.
- Do not twist or compress the power cord.
- Prior to usage, check the working voltage. In areas with unstable voltage, consider using a voltage stabilizer suitable for motor loads to ensure a stable voltage within the normal input range required by the machine:
110V/220V $\pm 10\%$.
- The device must be grounded properly to avoid electrical shock hazards. The device uses a three-pin plug, and the third pin is a ground pin. Only use earth/grounded sockets with a protective earth (PE) conductor.

8. Instructions for the Rotor and Tubes


- Insert the rotor to the motor shaft. Ensure the rotor is in position and connected to the shaft. Tighten the locking nut firmly by turning it clockwise with the wrench. Ensure the rotor is firmly secured on to the motor shaft without any slackness or movement. You can remove the rotor by turning the locking nut counterclockwise using the wrench.

- Attach and tighten the rotor lid.
- Insert tubes opposite each other in pairs into the rotor bores. To ensure symmetric loading, tubes that are arranged opposite each other must be of the same type and contain the same sample quantity. To minimize weight differences between filled sample tubes, we recommend balancing. This will reduce wear on the drive and also reduce operating noise.

9.Operation Settings

Operation Panel

①	Open/close the lid
②	Start/Stop program
③	Check the operation graph
④	Check the operation logs
⑤	Click to the settings menu
⑥	Back the the home screen
⑦	Set centrifugation time
⑧	Display the remaining centrifugation time
⑨	Set the speed
⑩	Display the real-time speed
⑪	Display the real-time temperature

Speed setting

Speed Setting

-

5000_{rpm}

+

30015000

Back

Save

Set the speed of centrifugation and adjust it by pressing “+/-” icons (increment is 100rpm) or use slide function for faster adjustment. The speed range is 300-15000rpm.

Time Setting

Time Setting

H	M	S	Range
98	08	58	30s-99h59m59s
99	09	59	
00	10	00	
01	11	01	
02	12	02	

Back

Save

Slide up and down to set the centrifugation time.


System Settings

System

Conversion

Device info

System settings

Brightness 50  100

Ramp rate

Up(1-3)  3  Down(1-3)  3 

Language ☐ 简体中文 ☒ English

Home

Settings

Logs

Graph

System

Conversion

Device info

System settings

RPM/RCF ☒ Speed (rpm)
☐ RCF(xg)

Temp. unit ☒ °C
☐ °F

Home

Settings

Logs

Graph

- Brightness: slide left or right to adjust the brightness.
- Continuous cooling: switch on and close the lid to start continuous cooling automatically.
- Three adjustable speed levels: speed can be accelerated from 0 to 15000rpm or decelerated from 15000rpm to 0 at three speed levels. When setting 1st, 2nd, 3rd speed level, the acceleration and deceleration time is 60s, 40s, 15s respectively.
- Language: Chinese or English
- RPM/RCF: speed or RCF
- Temp. Unit: °C or °F

Program Settings

Delete

Logs

No.	Status	Time	Description
1		2024-01-30 10:54:00	The device is on.

Home

Settings

Logs

Graph

New: Click “New” button to create a new program

Select: Choose a program to run (e.g. No.2) , then click “Select” button

Delete: Choose a program to delete (e.g. No.2) ,then click

“Delete” button

Edit: Select a program to edit (e.g. No.2), then click “Edit” button

10.Troubleshooting

Error messages


Problem	Cause	Solution
No display	No power connection	Check the power connection
	Damaged power socket	Replace the power socket
	Damaged power switch	Replace the power switch
	Others	Contact distributor or supplier
Centrifuge lid can't be opened.	The rotor is still running.	Wait for the rotor to stop.
	Power outage	1. Check the power supply 2. Activate the emergency lid release.
Loud noise during operation	The rotor is not loaded properly.	Switch off the device and tighten the rotor nut and rotor lid
	The tubes are asymmetrically loaded.	Switch off the device and load the tubes symmetrically.
Temperature display doesn't match the actual temperature.	Damaged temperature sensor.	Contact distributor or manufacturer
Shows “High Rotor Temperature. Device stopped automatically.”	The ambient temperature is too high.	The ambient temperature should be less than 35°C.

Note: In case of any error messages are not shown in above list, please contact your local distributor or manufacturer.

Emergency lid release

During power outage, the centrifuge lid cannot be opened using the “Open lid” button on the operation panel. In that case, the lid can be opened by the emergency lid release.

- Ensure the rotor is not running through the view port.
- Disconnect the device from the power supply.
- Insert the emergency lid release pin or a thin screwdriver into the emergency lid release hole and push it, then the lid will open automatically.

11. Maintenance and Cleaning

- Switch off the device and disconnect it from the power supply before starting cleaning or disinfection.
- Do not allow any liquids to penetrate the inside of the housing.
- When the centrifuge is not in use for an extended period, remove the rotor and apply rust inhibitor onto the rotor base (or drive shaft) to prevent rust. This ensures proper coordination

between the rotor base and rotor for balanced operation. To prevent rust, consider placing desiccants inside the centrifuge chamber to absorb moisture when needed.

- During use, the rotor may be contaminated due to damaged centrifuge tubes, overfilling, or improper sealing. After use, remove the rotor from the shaft, clean it promptly, and apply a protective agent. To prevent corrosion and associated risks, monitor sample volume and seal integrity, regularly inspect the seal ring for damage, and replace it promptly if aging or broken occurs to prevent leaks.
- Remove the rotor and clean it using a soft cloth soaked with diluted alcohol.
- Clean any dust deposits using a brush, but not a wire brush.
- Thoroughly clean the rotor and rotor chamber, and check the rotor for wear, damage and corrosion.
- Only use cleaning agents and disinfectants with pH between 6 and 8.
- If the device has been contaminated by aggressive chemicals, clean immediately using a mild cleaning agent.
- Do not exceed temperature of 121°C or a time of more than 20min. while autoclaving.
- Replace the seals of rotors after 50 autoclaving cycles. (Seal PN: 3.4.603440007)

12. Warranty

We guarantee that our scientific instruments adhere to the most rigorous engineering and quality standards. This instrument is warranted to be free from defects in materials and workmanship under normal use and service, for a period of 24 months from the date of dispatch. The warranty is extended only to the original purchaser. For claims under the warranty, please contact your local supplier. After the warranty period expires, the manufacturer retains the right to invoice the cost price for the repair or maintenance of a faulty device, along with any associated service fees.

Scope of Warranty

The following conditions are not covered under the warranty.

- Faults or damage caused by negligence, improper installation, improper operation, or failure to use and maintain the machine in accordance with the instructions in this operating manual.
- Issues caused by unauthorized disassembly or modification.

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