


Bird Gender Determination in 30min.

WizDx™ CrystalMix

Bird Gender

- High Sensitivity & Specificity
- Room temperature Storage/Shipping
- Ready-to-use Lyophilized PCR Kit
- Easy & Fast protocol



High Sensitive

The LoD of WizDx™ CrystalMix Bird Gender is 1.25 copies/ μ L.

Perfect Reproducibility

WizDx™ CrystalMix Bird Gender kits are precisely dispensed by fully automated dispensers in facilities certified with ISO 13485 and GMP. They are supplied in lyophilized form with a single test quantity, so all kits deliver exactly the same results.

Prevent Carryover-contamination

WizDx™ CrystalMix Bird Gender kit contains Uracil-DNA glycosylase (UDG) and dUTPs eliminating possible carryover contamination.

High Specificity

Accurate inspection is possible through the specific probe design for the Glycoprotein D gene region of HSV 1 & 2.

Easy & Accurate

Lyophilized PCR kits are very convenient to use because you can test immediately with the addition of sample DNA only, and provide more accurate results because you can drastically reduce experimental errors or contamination that may occur during the preparation process.

Convenience & Stable

Wizbiosolution's patented HNC(Heteromorphic Nanopolymer Coating) technology applied for WizDx™ CrystalMix PCR Kit enables the reagents stored and transported at room temperature with the benefits of user-convenience and stabilization.



Key Features of CrystalMix Bird Gender



Anti-Carryover Technology

Prevents cross-contamination



Hot-Start PCR Technology

Minimizes non-specific PCR reactions and increases specificity



Enhanced-stability technology

Sustains long-term storage at room temperature



Automated dispensing system

Maximizes reproducibility and minimizes experimental errors




Ready-to-use, Freeze-dried form

Allows convenient storage and minimal experimental process

Easy & Fast Molecular Diagnosis



Experimental Data



Ordering Information

Cat No.	Product	Package Size
DX0051	WizDx™ CrystalMix Bird Gender	96 T
CL0016	CLEO™ Q16 Real-time PCR system	1 set

Available sample, please contact to sales@wizbiosolution.com

